


Guide to Invasive Monitoring Configurations for Metavision



P1 (P1/P3)

Port 1 (Splitter Port Labeled As IP1) – Arterial

Port 3 (Splitter Port Labeled As IP2) – CVP (Central Venous Pressure – Mean Blood Pressure Lead 3)

P2 (P2/P4)

Port 2 (Splitter Port Labeled As IP1) – PA (Pulmonary Artery – Diastolic Blood Pressure Lead 2)

Port 4 (Splitter Port Labeled As IP2) – ICP (Intracranial Pressure – Mean Blood Pressure Lead 4)

